



Video-Med

By Juancarlos



My Goals


- 1- To Improve patients care outcomes
- 2- Reduce Cost
- 3- Educate Patient and provide experience
- 4- Deliver much more than a video conferencing
- 5- Delivering remote video consulting
- 6- Integrated medical devices
- 7- Easy to use for all patients

The first step in the process is to identify the problem. This involves gathering information about the situation and the people involved. Once the problem is identified, the next step is to analyze it. This involves breaking the problem down into its components and understanding how they are related. The third step is to develop a plan. This involves deciding on the best way to solve the problem and the resources needed to do so. The fourth step is to implement the plan. This involves putting the plan into action and monitoring progress. The final step is to evaluate the results. This involves assessing the effectiveness of the solution and making any necessary adjustments.



My purpose is to share via remote patient monitoring in one system with clinician. The device will include automatic video call routing by service with EKG and EMG sensors to process the electrical activity of the heart.

Visualize Concept




This is my original concept a screen in the center of the device with all the input sensors on the outside

Form



This was my form and
idea behind of the
video-med

Dimensions



All the dimensions are in inches. I started with the height 6", depth 2" width 4.5", screen width 3" and a height 2.5"



Video-Med

- Consist of the followings:
- Capacitive Touch Sensors as input
- An Actuator as output
- An Arduino
- An Olimexino 328
- Communication Network Protocol

Materials



Shield EKG



Shield EMG



Electrode Monitoring
medi face foam




3.2 TFT LCD Touch Screen



Arduino uno



Olimex Board




Prototype




--	--




LCD and Arduino Schematic



Breadboard + arduino + Olimexino328 and a 3.2 LCD screen



Schematic for Olimex with Arduino



Code for Olimex



ShieldEkgEmgCode.ino

```
#include <compat/deprecated.h>
#include <FlexiTimer2.h>
//http://www.arduino.cc/playground/Main/FlexiTimer2

// All definitions
#define NUMCHANNELS 6
#define HEADERLEN 4
#define PACKETLEN (NUMCHANNELS * 2 + HEADERLEN + 1)
#define SAMPFREQ 256 // ADC sampling rate 256
#define TIMER2VAL (1024/(SAMPFREQ)) // Set 256Hz sampling frequency
#define LED1 13
#define CAL_SIG 9

// Global constants and variables
volatile unsigned char TXBuf[PACKETLEN]; //The transmission packet
volatile unsigned char TXIndex; //Next byte to write in the transmission packet.
volatile unsigned char CurrentCh; //Current channel being sampled.
volatile unsigned char counter = 0; //Additional divider used to generate CAL_SIG
volatile unsigned int ADC_Value = 0; //ADC current value

//~~~~~
// Functions
//~~~~~

/*****
 * Function name: Toggle_LED1 */
 * Parameters */
 * Input : No */
 * Output : No */
 * Action: Switches-over LED1. */
*****/
void Toggle_LED1(void){

if((digitalRead(LED1))==HIGH){ digitalWrite(LED1,LOW); }
else{ digitalWrite(LED1,HIGH); }

}
```

```

    void toggle_GAL_SIG(void){
        #if(digitalRead(CAL_SIG) == HIGH){ digitalWrite(CAL_SIG, LOW);}
        else{ digitalWrite(CAL_SIG, HIGH);}
    }

}

/*-----*/
/* Function name: setup */
/* Parameters */
/* Input : No */
/* Output : No */
/* Action: Initializes all peripherals */
/*-----*/
void setup() {

    noInterrupts(); // Disable all interrupts before initialization

    // LED1
    pinMode(LED1, OUTPUT); //Setup LED1 direction
    digitalWrite(LED1,LOW); //Setup LED1 state
    pinMode(CAL_SIG, OUTPUT);

    //Write packet header and footer
    TXBuf[0] = 0xa3; //Sync 0
    TXBuf[1] = 0x5a; //Sync 1
    TXBuf[2] = 2; //Protocol version
    TXBuf[3] = 0; //Packet counter
    TXBuf[4] = 0xc2; //CH1 High Byte
    TXBuf[5] = 0x00; //CH1 Low Byte
    TXBuf[6] = 0x02; //CH2 High Byte
    TXBuf[7] = 0x00; //CH2 Low Byte
    TXBuf[8] = 0x02; //CH3 High Byte
    TXBuf[9] = 0x00; //CH3 Low Byte
    TXBuf[10] = 0x02; //CH4 High Byte
    TXBuf[11] = 0x00; //CH4 Low Byte
    TXBuf[12] = 0x02; //CH5 High Byte
    TXBuf[13] = 0x00; //CH5 Low Byte
    TXBuf[14] = 0x02; //CH6 High Byte
    TXBuf[15] = 0x00; //CH6 Low Byte
    TXBuf[2 * NUMCHANNELS + HEADERLEN] = 0x01; // Switches state

    // Timer2
    // Timer2 is used to setup the analog channels sampling frequency and packet update.
    // Whenever interrupt occurs, the current read packet is sent to the PC
    // In addition the CAL_SIG is generated as well, so Timer1 is not required in this case!
    FlexTimer2::set(TIMER2VAL, Timer2_Overflow_ISR);
    FlexTimer2::start();

    // Serial Port
    Serial.begin(57600);
    //Set speed to 57600 bps

    // MCU sleep mode = idle.
    //outb(MCUCR,(inp(MCUCR) | (1<<SE) & ~(1<<SM0) | ~(1<<SM1) | ~(1<<SM2)));

    interrupts(); // Enable all interrupts after initialization has been completed
}

/*-----*/
/* Function name: Timer2_Overflow_ISR */
/* Parameters */
/* Input : No */
/* Output : No */
/* Action: Determines ADC sampling frequency. */
/*-----*/
void Timer2_Overflow_ISR()
{
    // Toggle LED1 with ADC sampling frequency /2
    Toggle_LED1();


    //Read the 6 ADC inputs and store current values in Packet
    for(CurrentCh=0;CurrentCh<6;CurrentCh++){
        ADC_Value = analogRead(CurrentCh);
        TXBuf[((2*CurrentCh) + HEADERLEN)] = ((unsigned char)(ADC_Value & 0xFF00) >> 8); // Write High Byte
        TXBuf[((2*CurrentCh) + HEADERLEN + 1)] = ((unsigned char)(ADC_Value & 0x00FF)); // Write Low Byte
    }

    // Send Packet
    for(TXIndex=0;TXIndex<17;TXIndex++){
        Serial.write(TXBuf[TXIndex]);
    }

    // Increment the packet counter
    TXBuf[3]++;

    // Generate the CAL_SIGNAL // increment the divider counter
    counter++;
    if(counter == 12){ // 250/12/2 = 10.4Hz ->Toggle frequency
        counter = 0;
        toggle_GAL_SIG(); // Generate CAL signal with frequ ~10Hz
    }
}

```

```
■ /******  
■ /* Function name: loop */  
■ /* Parameters */  
■ /* Input : No */  
■ /* Output : No */  
■ /* Action: Puts MCU into sleep mode. */  
■ /******  
■ void loop() {  
■  
■ __asm__ __volatile__ ("sleep");  
■  
■ }
```

Code for Processing



ekg_processing_code.txt

- import processing.serial.*;
- Serial myPort; // The serial port
- int counter = 0;
- int inByte = 0;
- float x1, y1, x2, y2, y;
- float spacing = 1.0;
- color RED = color(255, 0, 0);
- color GREEN = color(0, 255, 0);
- color BLUE = color(12, 16, 255);
- color WHITE = color(255, 255, 255);
- color BLACK = color(0, 0, 0);
- void setup () {
- size(800, 400); // window size
- // **** List available serial ports **** //
- String[] ports = Serial.list();
- print("No. ports = ");
- println(ports.length);
- for (int i = 0; i < ports.length; i++) {
- print "[" + i + "];
- println(ports[i]);
- }
- // **** Enter selected port number here **** //
- String portName = Serial.list()[0];
- print("port selected = ");
- println(portName);
- myPort = new Serial(this, portName, 9600);
- background(BLUE);
- }

- void serialEvent (Serial myPort) {
- int inByte = myPort.read();
- println(inByte);
- y = height - inByte;
- }


- void draw () {
- // line color:
- stroke(WHITE);
-
- if (counter > width/spacing) {
- counter = 0;
- background(BLUE);
- }
-
- if (counter == 0) {
- x1 = 0;
- y1 = y;
- }
-
- if (counter > 0) {
- x2 = counter*spacing;
- y2 = y;
- line(x1, y1, x2, y2);
- x1 = x2;
- y1 = y2;
- }
- counter++;
- }



Instructable

■ INTRODUCTION:

- Video med uses information Technology and telecommunication to provide clinical health care from a distance. This Technology allows communication between patients and medical personnel and the transmission of health informatics from one site to another.



1- Materials

- [Arduino Uno](#)



- [ECG-PRO-3-WAY-CABLE](#)




- [Olimex Shield Ecg, Ekg](#)




- [SHIELD-EKG-EMG-PA](#)



- [Breadboard](#)



- [3.2 sainsmart tft lcd](#)



- [ECG-GEL-ELECTRODE](#)



- [TFT Shield For Arduino UNO R3](#)



3D Models



telemedprototype.123dx




telemed print front hole.stl



telemed 1 print.stl

Schematics

LCD and Arduino Schematic




olimexarduino.sch

Arduino Code



sketch_dec06jca.ino

Processing Code



ekg_processing_ino.txt

Step by Step Desing

- 1) My first step was to prepare the [Olimex](#) 328 Board by removing its wire between the R6 pads. Failure to perform this step could result in electrical failure from increase power consumption.



- 2) I downloaded the IDE package from the [Arduino](#) website.
- 3) I downloaded 2 libraries required to run Olimex ([Timer one](#), [Timer2](#))
- 4) Next I arranged the EKG shield jumpers as follow:
 - REF_E – closed
 - 3.3V/5V – 5V position
 - D4/D9 – D9 position
 - ANI_SEL – 1 position

- 5) I connected the Shield to The Arduino Board.




- 6) I then connected the Arduino Olimex board to the USB
- 7) I then Downloaded the Code Shield EkgrmgCode. I compiled the code and I rang it.
- 8) I then downloaded the code for EKG Processing

- 9) I Then attached the Ekg gel pads to “the patient”



- 10) The Ekg Processing code was compile and rang.
- 11) The screen Opened and the Electro Cardiogram (EKG) was displayed.



Difficulties encountered

- 1) In setting up my initial LCD screen, I learned that the 3.2 LCD screen was compatible with Arduino mega and not Arduino Uno
- 2) My LCD screen no longer supported myGLCD library. I then researched and found an alternate library [Open.myGLCD](#)
- 3) Manipulating the code for showing the EKG on the LCD screen did not work. This need more time to solve with additional trial and error.

What I learned

What I learned is how to manipulate LCD screens, how to trouble shoot, including testing and debugging. The process was fascinating, challenging but very rewarding. It was valuable experience to go from concept inspiration (Ideation) to prototyping to 3D Modeling, Designing schematics and finally coding.

What is still Outstanding

While I was able to show a graph on the LCD screen, I was unable to debug the code to show the actual EKG on the LCD Screen.